

PANALYTICAL X'PERT PRO: Quick Start Guide

The X'Pert Pro is a versatile X-ray Diffraction (XRD) system with interchangeable "PreFIX" modules. The default configuration enables basic measurements, such as θ - 2θ scans for phase identification. The standard training for the X'Pert Pro does not authorize you to change the PreFIX modules or the Sample Stage. If you require a different configuration, contact ADL staff to make the change for you.

1. Starting up the system

1. Log onto the instrument PC.
2. Open the **Data Collector** software.
3. On the top toolbar, select "**Instrument**", then select "**Connect**" in the drop-down.
 - a. On the pop-up, select the "**USER**" configuration and "**Triple Axis / Coll. 0.18**" radio button, then click **OK**.
 - b. Another pop-up will alert you to the assumed instrument settings. Click **OK**.
4. Once connected, check the toolbar that has appeared on the left.
 - a. On the "**Incident Beam Optics**" tab, set the Beam attenuator usage to "**Do not switch**" and select the "**Activated**" checkbox (if it is not already).
 - b. On the "**Instrument Settings**" tab, set the Generator to **45 kV / 40 mA**.

If the XRD hasn't been used for more than 24 hours, you may want to run the "fast" (~6 min) X-ray gun breeding procedure via the button in the "X-ray Settings" pop-up.

2. Mounting samples

Solid samples can be mounted using double-sided tape or, if they are light and thin enough, the three sliding clips on the basic sample stage.

Powder samples must be consolidated in some way; either by spreading a uniform layer onto double-sided tape, pressing into a solid pellet using a pellet press, etc.

You may want to mount your sample on an amorphous material, such as a glass slide, to reduce potential spectral contamination from the aluminum mounting surface.

3. Aligning samples

1. Ensure the Beam attenuator usage is set to "**Do not switch**" and "**Activated**".
 2. Ensure the sample stage is at its home position (all stage positions set to zero).
 3. Ensure the beam optics are appropriate for alignment.
 - a. Insert the appropriately sized **vertical beam mask** into the incident beam optics.
 - b. Insert the **1/32° divergence slit** into the incident beam optics.
 - c. Insert the **0.18° parallel plate collimator slit** into the diffracted beam optics.
-
4. Begin by aligning **2θ**.
 - a. On the top toolbar, select "**Measure**", then select "**Manual Scan**" in the drop-down.
 - i. On the pop-up, select "**2Theta**" as the Scan Axis.
 - ii. Type in a Range of **1°**, Step size of **0.01°**, and Time per step of **0.2 s**.
 - iii. Then select "**Start**" and wait for the scan to finish.
 - b. Right click on the completed scan and select "**Move Mode**" from the context menu.
 - c. Click and drag the **vertical green line** onto the peak maximum.
 - d. Close the scan.
-
5. Now, manually move the **Z** position to approximately the correct position for your sample thickness (typically somewhere between 6 - 10 mm).
 6. Next, align **Z**.
 - a. On the top toolbar, select "**Measure**", then select "**Manual Scan**" in the drop-down.
 - i. On the pop-up, select "**Z**" as the Scan Axis.
 - ii. Type in an appropriate Range, Step size, and Time per step. Try to pick values that keep the total scan time under 1 minute.
 - iii. Then select "**Start**" and wait for the scan to finish.
 - b. Right click on the completed scan and select "**Move Mode**" from the context menu.
 - c. Click and drag the **vertical green line** onto the midpoint of the intensity dropoff.
 - d. Close the scan.
-

7. Finally, align Ω .
 - a. On the top toolbar, select "**Measure**", then select "**Manual Scan**" in the drop-down.
 - i. On the pop-up, select "**Omega**" as the Scan Axis.
 - ii. Type in an appropriate Range, Step size, and Time per step. Again, try to keep the total scan time under 1 minute. Typically you will need a Range greater than 3 - 4°.
 - iii. Then select "**Start**" and wait for the scan to finish.
 - b. Right click on the completed scan and select "**Move Mode**" from the context menu.
 - c. Click and drag the **vertical green line** onto the peak/midpoint.
 - d. Close the scan.
-

8. Repeat the **Z** and Ω alignment procedures until you are satisfied with the alignment.
 9. On the top toolbar, select "**User Settings**", then "**Sample Offsets**" in the drop-down.
 - a. Set the Current Position for both **2 θ** and Ω to be zero.
 - b. Click "**OK**".
-

Your sample should now be properly aligned.

4. Taking data

1. Start by ensuring the system settings are appropriate for taking data.
 - a. On the "**Incident Beam Optics**" tab, set the Beam attenuator usage to "**At preset intensity**".
 - b. Insert the desired **divergence slit** into the incident beam optics; e.g. 1/4°, 1/2°, etc.
 - c. If you don't need the best possible resolution, it is recommended to remove the **0.18° parallel plate collimator** slit to improve signal intensity.
2. Next, generate a new program that defines your **θ -2 θ** scan.
 - a. On the top toolbar, select "**File**", then select "**New Program**" in the drop-down.
 - b. On the pop-up, select "**Absolute Scan**" in the Program type drop-down.
 - c. Click "**OK**".

- d. On the next window, select "**2Theta-Omega**" in the Scan Axis drop-down.
 - e. Enter your desired scan parameters; i.e. start/end angles, step size, time per step.
 - f. Close the window using the red **X** button in the upper right corner. The **Data Collector** software will then ask if you want to save the program you just created. Save it in your personal program folder.
3. Finally, run your saved program.
- a. On the top toolbar, select "**Measure**", then select "**Program**" in the drop-down.
 - i. Click the "**Browse**" button.
 - ii. In the file explorer, select the file containing your previously saved program.
 - b. On the next window, click the button with an icon that resembles a file folder.
 - i. In the file explorer, select the location you want to save the finished data file.
 - ii. Type in a filename, then click "**OK**".
 - c. Type in any additional metadata you want saved to the file; i.e. ID, name, etc.
 - d. Click "**OK**".

Your scan should now be running! Wait for it to finish.

5. Shutting down the system

1. Return the system to its standby settings.
 - a. On the "**Incident Beam Optics**" tab, set the Beam attenuator to "**Do not switch**" and select the "**Activated**" checkbox.
 - b. On the "**Instrument Settings**" tab, set the Generator to **40 kV / 20 mA**.
2. Remove your previously set offsets.
 - a. On the top toolbar, select "**User Settings**", then "**Sample Offsets**" in the drop-down.
 - i. Click "**Clear All Offsets**" and then click "**OK**".
3. Return the sample stage back to its home position (all motors at zero).
4. On the top toolbar, select "**Instrument**", then select "**Disconnect**" in the drop-down.
5. You can now close the **Data Collector** software.